

OneRoof: A Pipeline Prototype for Base-, Variant-, and Consensus-calling under One Proverbial Roof

Table of contents

Overview	1
Quick Start	2
Detailed Setup Instructions	3
Configuration	3
Developer Setup	5
Contributing	6
Is it any good?	6
Citation	6

Overview

oneroof is a pipeline designed to take a common series of bioinformatic tasks (see below) and put them under “one roof”. We mean this quite literally: the pipeline will perform at its best when run on networked devices in the same building.

oneroof was originally developed in the early stages of the United States Bovine Highly Pathogenic Avian Influenza (HPAI) outbreak of 2024, when we wanted one, configurable, easy-to-run pipeline that would do all of the following seamlessly:

1. Handle super-accuracy basecalling with GPU acceleration on pod5-formatted Nanopore signal files, working on GCP, AWS, or Slurm if need be.
2. Demultiplex BAM-formatted reads that come out of basecalling.
3. Perform the above two steps as signal files become available, either locally or remotely via an SSH client.
4. Accept raw read BAMs or FASTQs if basecalling and demultiplexing have already been performed elsewhere.
5. Accept paired Illumina reads in addition to Nanopore reads.
6. Use forward and reverse primer sequences to select only those reads that represent complete amplicons.
7. Trim away primers—and also any bases that are upstream of the forward primer or downstream of the reverse primer, *while only retaining those reads identified as complete amplicons by containing both primers*.
8. Align to a custom reference with the proper presets for the provided data.
9. Call variants and consensus sequences with appropriate settings for the provided data.
10. Perform tree building with `nextclade`, quality introspection with `multiQC`, and error correction based on the input sequence platform.

Though many excellent pipelines currently exist, e.g. `nf-core/viralrecon`, `epi2me-labs/wf-amplicon`, and `nf-core/nanoseq`, none of these pipelines quite handled all of the above. `oneroof` seeks to handle these requirements while remaining highly configurable for users, highly modular for developers, and easy to control in the command line for both.

Overall, `oneroof` can be summarized as a variant-calling pipeline written in and managed by Nextflow. Its software dependencies are provided through containers or through an environment assembled by `pixi`. To run it on your own Nanopore pod5s with Docker containers, simply run something like:

```
nextflow run nrminor/oneroof \
--pod5_dir my_pod5_dir \
--primer_bed my_primers.bed \
--refseq my_ref.fasta \
--ref_gbk my_ref.gbk \
--kit "SQK-NBD114-24"
```

These are the core elements required to run on Nanopore data: a directory of pod5 files, a BED file of primer coordinates, a reference sequence in FASTA and Genbank format, and the Nanopore barcoding kit used.

And for Illumina paired-end reads, it's even simpler:

```
nextflow run nrminor/oneroof \
--illumina_fastq_dir my_illumina_reads/
```

If you want to use Apptainer containers instead of Docker, just add `-profile apptainer` to either of the above `nextflow run` commands. And if you don't want to use containers at all, simply run `pixi shell --frozen` to bring all the pipeline's dependencies into scope and then add `-profile containerless` to your `nextflow run` command.

Nextflow pipelines like this one have a ton of configuration, which can be overwhelming for beginners and new users. To make this process easier, we're developing a Terminal User Interface (TUI) to guide you through setup. Please stay tuned!

Quick Start

For most users, `oneroof` will have two core requirements: The Docker container engine, available [here](#), and Nextflow, available [here](#). For users interested in super-accuracy basecalling Nanopore signal files, an on-board GPU supported by the Dorado basecaller is also required.

All remaining software dependencies will be supplied through the pipeline's Docker image, which will be pulled and used to launch containers automatically.

From there, the pipeline's three data dependencies are sequence data in BAM, FASTQ, or POD5 format, a BED file of primer coordinates, and a reference sequence in FASTA and Genbank for-

mat. For Nanopore data, a barcoding kit identifier is also required. Simply plug in these files to a command like the above and hit enter!

Detailed Setup Instructions

Configuration

Most users should configure oneroof through the command line via the following parameters:

Command line argument	Default value	Explanation
<code>--primer_bed</code>	None	A bed file of primer coordinates relative to the reference provided with the parameters <code>refseq</code> and <code>ref_gbk</code> .
<code>--refseq</code>	None	The reference sequence to be used for mapping in FASTA format.
<code>--ref_gbk</code>	None	The reference sequence to be used for variant annotation in Genbank format.
<code>--remote_pod5_location</code>	None	A remote location to use with a TCP stream to watch for pod5 files in realtime as they are generated by the sequencing instrument.
<code>--file_watcher_config</code>	None	Configuration file for remote file monitoring. An example can be found at <code>conf/file_watcher.template.yml</code> .
<code>--pod5_dir</code>	None	If a remote pod5 location isn't given, users may provide a local, on-device directory where pod5 files have been manually transferred.
<code>--precalled_staging</code>	None	A local directory to watch for Nanopore FASTQs or BAMs as they become available. This is for cases where basecalling is being performed via another workflow, resulting in BAMs or FASTQs that are gradually transferred into <code>params.precalled_staging</code> as basecalling completes.
<code>--prepped_data</code>	None	If pod5 files have already been basecalled, demultiplexed, and transferred to a local directory accessible to oneroof users can specify their location with <code>--prepped_data</code> .

Command line argument	Default value	Explanation
<code>--illumina_fastq_dir</code>	None	If users have Illumina data to be processed, they may specify their paired-end FASTQ files' location with <code>illumina_fastq_dir</code> .
<code>--model</code>	sup@latest	the Nanopore basecalling model to apply to the provided pod5 data (defaults to the latest super-accuracy version)
<code>--kit</code>	None	The Nanopore barcoding kit used to prepare sequencing libraries.
<code>--pod5_batch_size</code>	all pod5s	How many pod5 files to basecall at once. With a single available GPU, all pod5 files should be basecalled together, so this parameter defaults to telling Nextflow to take all pod5 files at once.
<code>--basecall_max</code>	1	If basecalling pod5 files is to be parallelized across multiple available GPUs, this parameter tells Nextflow how many parallel instances of the basecaller to run at once (defaults to 1).
<code>--max_len</code>	12345678910	The maximum acceptable length for a given read.
<code>--min_len</code>	1	The minimum acceptable length for a given read.
<code>--min_qual</code>	0	The minimum acceptable average quality for a given read.
<code>--max_mismatch</code>	0	The maximum number of mismatches to allow when finding primers
<code>--downsample_to</code>	0	Desired coverage to downsample to, with a special value of 0 to indicate 0 downsampling
<code>--secondary</code>	None	Whether to turn on secondary alignments for each amplicon. Secondary alignments can increase depth at the cost of more reads potentially mapping to the wrong locations. By default, secondary alignments are off.
<code>--min_consensus_freq</code>	0.5	The minimum required frequency of a variant base to be included in a consensus sequence.

Command line argument	Default value	Explanation
<code>--min_depth_coverage</code>	20	Minimum required depth of coverage to call a variant.
<code>--results</code>	results/	Where to place results.
<code>--cleanup</code>	false	Whether to cleanup work directory after a successful run.

Note that oneroof checks for how to gather data in a particular order for Nanopore data, which mirrors the order they are listed in the table above:

1. First, oneroof checks if the user has supplied a remote location with `--remote_pod5_location`, at which point it will launch the file watcher module and begin transferring and basecalling batches of pod5 files. The size of these batches can be controlled with `--pod5_batch_size`, which defaults to basecalling batches of 100 pod5s.
2. If no remote directory is provided, oneroof checks for a local pod5 directory from `--pod5_dir`.
3. If the user doesn't provide a local `--pod5_dir`, oneroof assumes that pre-basecalled BAMs or FASTQs are being provided. These can either be watched for with the directory from `--precalled_staging`, or run immediately with the directory from `--prepped_data`.

Developer Setup

oneroof depends on software packages supplied through various conda registries as well as through PyPI, the Python Package Index. To unify these various channels, we used the relatively new `pixi` package and environment manager, which stores all dependencies from both locations in the file `pyproject.toml`.

To reproduce the environment required by this pipeline, make sure you are on a Mac, a linux machine, or a Windows machine using Windows Subsystem for Linux. Then, to reproduce the environment, install `pixi` with:

```
PIXI_ARCH=x86_64 curl -fsSL https://pixi.sh/install.sh | bash
```

Download the pipeline with:

```
git clone https://github.com/nrminor/oneroof.git && cd oneroof
```

And then open a `pixi` subshell within your terminal with:

```
pixi shell --frozen
```

As long as you are using a supported system, the pipeline should run within that subshell. You can also run the pipeline within that subshell without containers using the “containerless” profile:

```
nextflow run . \
-profile containerless \
--pod5_dir my_pod5_dir \
--primer_bed my_primers.bed \
--refseq my_ref.fasta \
--ref_gbk my_ref.gbk \
--kit "SQK-NBD114-24"
```

Especially on Apple Silicon Macs, this will reduce the overhead of using the Docker Virtual Machine and allow the pipeline to invoke tools installed directly within the local project environment.

Note also that more information on the repo's files is available in our developer guide.

Contributing

Contributions, feature requests, improvement suggestions, and bug reports via GitHub issues are all welcome! For more information on how to work with the project and what all the repo files are, see our developer guide.

Is it any good?

Yes.

Citation

Coming soon!